



ROTARY SEALER

[Code: SSAUM086]



SEALING MACHINE

Completion of a Sterile Barrier System.
Extremely easy to use. Solid, durable.
Polished chromium frame, hygienic,
No confusing buttons

Once power is on and temperature set, it works when it is needed and stops automatically when the work is done.
Accurate temperature control allows to seal Tyvek pouches

- Corresponds to DIN 58953-7
- Constructed from Stainless Steel for hygienic use.
- Multi line seal
- Temperature control with a microprocessor.
- Temperature tolerance, ± 1
- Automatic in-feed with optic sensor.
- 4 Language
- Temperature can be adjusted between 50°C to 200°C.
A wide variety of packaging materials can be hermetically sealed.
- Menu driven device through digital LCD display
- 8-12 Meter/minute speed
- Adjustable seal distance
- Automatic switch-off when temperature deviates



ROLL DISPENSER&CUTTER

CODE	DIMENSIONS
SS2001	500 x 620 x 500 mm
SS2002	500 x 102 x 250 mm

- Easy to use
- Can be mounted on a wall or can also be used on a work table, desk...
- Made of Stainless Steel
- Strong durable
- Cutting system do not generate paper dust.
- Almost lifetime rotary knife.
- Knife sharpens automatically while working

