

WITH  
PRINTING UNIT

# ROTARY SEALER

with Printing Unit



- Large LCD screen.
- Built in printer: Seal date, Expiry date, Operator name/code,
- Batch Counter
- Elapsed Hour Meter
- Clock and Calendar
- Data storage, when the device is switched off
- Automatically updated current and expiry date
- Globally understandable coding. No need to select language
- Adjustable font width
- Printer can be set off and on
- Optic sensor for start-up
- Sleep mode when not used.
- Corresponds to DIN 58953-7
- Constructed from Stainless Steel for hygienic use.
- Multi line seal
- Temperature control with a microprocessor.
- Temperature tolerance,  $\pm 1\%$
- Temperature can be adjusted between 50°C to 200°C. A wide variety of packaging materials can be hermetically sealed.
- 8-12 Meter/minute speed
- Adjustable seal distance
- Automatic switch-off when temperature deviates